


NEW PRODUCT



Katie Chair
not included

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

KATIE

SIDE FOLDING TABLE

RG63SFT

FEATURES

Side folding table can be fitting to the:

- Katie High Back Chair (RG635)
- Katie Low Back Chair (RG630)

Can be fitted on left or right side of Katie Chair
Ideal for holding small items
(Drink bottle / Phone / Book / TV remote etc)
Hard wearing laminated platform with plastic edging
400 x 200mm table platform
Table can be folded down when not in use
Easy to engage / disengage folding table
Large user friendly release lever
Strong lightweight construction

SPECIFICATIONS (All dimensions shown are nominal)

Dimensions:

Table Width: **200mm**

Table Depth: **400mm**

Maximum Weight Loading: **10kg** (evenly distributed)

Product Weight: **3kg**



Please contact **REDGUM Brand** if you require further information on this product.

P: 08 9248 4180

E: orders@for-de.com.au

REDGUM brand - 1 Business way, Malaga, Western Australia 6090
Part of the FOR-DE Group Pty Ltd

redgumbrand.com.au