SHOPRIDER mobility

LION TILT TXR





FEATURES

G424TXR

- Mid-wheel drive for maximum maneuverability
- 6-wheel design for ultimate stability
- Full featured suspension system (Elastomer suspension)
- Programmable joystick controller - VR2
- Deluxe, heavy duty reclining captain seat
- Adjustable footplate and arm assemblies
- Flat-free, non-marking tyres
- 8A off-board battery charger
- Power tilt seating (REAC)
- High ground clearance with 14" drive wheel
- Certified safety belt with rated bolt attachment

SPECIFICATIONS

- Base Dimensions : **1069mm (L) x 625mm (W)**
- Standard Seat: **18" x 18"**
- Ground Clearance: **69.9mm**
- Obstacle Climbing: **40mm**
- Maximum range: **Up to 25km**
- Safe climb angle: **6°**
- Batteries: **2 x 12v50Ah (supplied)**
- Battery Weight: **17.7kg**
- Charger: **8 Amp**
- Product Weight: **98.4kg - without seat & batteries**
- Maximum User Weight: **136kg - on flat ground**

OPTIONS

18 X 18 PAN SEAT

20 X 20 PAN SEAT

18 X 20 PAN SEAT

POWER ELEVATED FOOTREST

NOTE: These pan seat options may be suitable for your specific seating solution - enquire with your local Authorised Dealer for a trial or to learn more.

BURGUNDY RED

Scan the QR code to visit the **Shoprider website**

Please contact **SHOPRIDER** for more information on this product. This product sheet is subject to change and **SHOPRIDER** can alter it at any time.

CRASH TESTED
ISO 7176:19

P: **08 9248 4180**
E: **orders@for-de.com.au**

ASK FOR THE BRAND

shoprider.com.au

IMPORTANT NOTE: It is the responsibility of the user to ensure the suitability of this product for its intended use.
The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

SHOPRIDER strives to provide the highest quality products possible that are subject to change without notice.